

Appendix (A): Clinic Information

Please fill in the blanks below and use this as a Fax cover

Clinic Name		
Lab director's Name		
Operator Name		
E-mail		
Clinic's address		
Clinic's number	Tel)	Fax)
FREND Serial number	Setup <input type="checkbox"/> System View <input type="checkbox"/> Serial Number F10U	
Check in the box of the item that you installed in your clinic.	Select installed items	Appendix included in this FAX
	<input type="checkbox"/> Free T4	
	<input type="checkbox"/> PSA	
	<input type="checkbox"/> Testosterone	
	<input type="checkbox"/> TSH	
	<input type="checkbox"/> 25-Hydroxy Vitamin D	

Contact information

Address	220 Bear Hill Road, STE 102. Waltham, MA 02451
E-mail	support@nanoentek.com
FAX	(781) 790-5649
Tel.	781-472-2558

Appendix (B): Verification of Reportable Range / Calibration

Materials to use:

1. FREN[™]D Linearity Kit from Nova-One Diagnostics (Nova-One.org)
2. FREN[™]D Cartridges (fT4 Reagents), FREN[™]D System and Auto Processor (AP)
3. Pipette and tips: 70 µL placed in Pretreatment Tube, place fT4 cartridge and Pretreatment tube in AP.
When completed, place cartridge from AP into FREN[™]D. Use LIN mode.

Assay Name	Cartridge LOT#	Exp. Date (MM/DD/YR)
fT4		/ /

Verification of Reportable Range

Nova-One Lot #: 6210A21002

Date of testing (MM/DD/YR): / /



Calibration Verification Results

(Unit: ng/dL)

	Vial - 1	Vial - 5	
Range	<0.40 – <0.40 – 1.21	2.47 – 4.11 – 5.75	
Run #1			
Run #2			
Run #3			
	Vial - 2	Vial - 3	Vial - 4
Range	1.26 – 2.10 – 2.94	1.85 – 3.09 – 4.33	2.20 – 3.66 – 5.12
Run #1			
Run #2			
Run #3			

* Run All Vials in Triplicate on same day
 * Repeat if the result is out of range.
 * Copy additional pages as needed in case for more items or if you have another FREN[™]D System.
 * Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (B): Verification of Reportable Range / Calibration

Materials to use:

1. FREND™ Linearity Kit from Nova-One Diagnostics (Nova-One.org)
2. FREND™ Cartridges (PSA Reagents), FREND™ System
3. Pipette and tips: 35 µL placed directly onto PSA Cartridge. Place PSA cartridge on FREND. Use LIN mode.

Assay Name	Cartridge LOT#	Exp. Date (MM/DD/YR)
PSA		/ /

Verification of Reportable Range

Nova-One Lot #: 6210A21002

Date of testing (MM/DD/YR): / /



Calibration Verification Results

(Unit: ng/mL)

	Vial - 1	Vial - 5	
Range	0.09 – 0.23 – 0.38	>25.00	
Run #1			
Run #2			
Run #3			
	Vial - 2	Vial - 3	Vial - 4
Range	3.84 – 6.40 – 8.96	7.06 – 11.77 – 16.48	13.04 – 21.73 →25.00
Run #1			
Run #2			
Run #3			

* Run All Vials in Triplicate on same day
 * Repeat if the result is out of range.
 * Copy additional pages as needed in case for more items or if you have another FREND™ System.
 * Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (B): Verification of Reportable Range / Calibration

Materials to use:

1. FREN[™]D Linearity Kit from Nova-One Diagnostics (Nova-One.org)
2. FREN[™]D Cartridges (Testosterone Reagents), FREN[™]D System and Auto Processor (AP)
3. Pipette and tips: 70 µL placed in Pretreatment Tube, place Testosterone cartridge and Pretreatment tube in AP. When completed, place cartridge from AP into FREN[™]D. Use LIN mode.

Assay Name	Cartridge LOT#	Exp. Date (MM/DD/YR)
Testosterone		/ /

Verification of Reportable Range

Nova-One Lot #: 6210A21002

Date of testing (MM/DD/YR): / /



Calibration Verification Results

(Unit: ng/dL)

	Vial - 1	Vial - 5	
Range	<20.00 – <20.00 – 33.02	683.3 – 1139 – >1500	
Run #1			
Run #2			
Run #3			
	Vial - 2	Vial - 3	Vial - 4
Range	136.7 – 227.8 – 318.9	253.9 – 423.1 – 592.3	468.3 – 780.5 – 1093
Run #1			
Run #2			
Run #3			

* Run All Vials in Triplicate on same day

* Repeat if the result is out of range.

* Copy additional pages as needed in case for more items or if you have another FREN[™]D System.

* Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (B): Verification of Reportable Range / Calibration

Materials to use:

1. FREN[™]D Linearity Kit from Nova-One Diagnostics (Nova-One.org)
2. FREN[™]D Cartridges (TSH Reagents), FREN[™]D System
3. Pipette and tips: 35 µL placed directly onto TSH Cartridge. Place TSH cartridge on FREN[™]D. Use LIN mode.

Assay Name	Cartridge LOT#	Exp. Date (MM/DD/YR)
TSH		/ /

Verification of Reportable Range

Nova-One Lot #: 6210A21002

Date of testing (MM/DD/YR): / /



Calibration Verification Results

(Unit: mIU/L)

	Vial - 1	Vial - 5	
Range	<0.06 – <0.06 – 0.26	22.75 – >25.00 – >25.00	
Run #1			
Run #2			
Run #3			
	Vial - 2	Vial - 3	Vial - 4
Range	3.36 – 5.59 – 7.83	6.09 – 10.14 – 14.20	10.09 – 16.81 – 23.54
Run #1			
Run #2			
Run #3			

- * Run All Vials in Triplicate on same day
- * Repeat if the result is out of range.
- * Copy additional pages as needed in case for more items or if you have another FREN[™]D System.
- * Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (B): Verification of Reportable Range / Calibration

Materials to use:

1. FREND™ Linearity Kit from Nova-One Diagnostics (Nova-One.org)
2. FREND™ Cartridges (Vitamin D Reagents), FREND™ System and Auto Processor (AP)
3. Pipette and tips: 35 µL placed into dilution buffer and mixed it. 70 µL from the mixture on the pretreatment tube, place the Vitamin D cartridge and pretreatment tube in AP. When complete, place the cartridge from AP into FREND. Use LIN mode

Assay Name	Cartridge LOT#	Exp. Date (MM/DD/YR)
Vitamin D		/ /

Verification of Reportable Range

Nova-One Lot #: 6210A21002

Date of testing (MM/DD/YR): / /



Calibration Verification Results

(Unit: ng/mL)

	Vial - 1	Vial - 5	
Range	<13.00 – <13.00 – 23.24	93.80 – >96.00 – >96.00	
Run #1			
Run #2			
Run #3			
	Vial - 2	Vial - 3	Vial - 4
Range	<13.00 – 26.35 – 42.30	27.30 – 45.50 – 63.70	54.02 – 90.04 – >96.00
Run #1			
Run #2			
Run #3			

* Run All Vials in Triplicate on same day

* Repeat if the result is out of range.

* Copy additional pages as needed in case for more items or if you have another FREND™ System.

* Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (C): Verification of Accuracy / Precision

Free T4

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean	Expected Ranges
Level 1	0.97	0.53 – 1.40
Level 2	1.93	1.02 – 2.84
LOT #	6361A23002, 6362A23002	
Exp. Date	2025-11-30	

* Specifications vary by Lot#. For lot-specific values, please visit NovaOne’s website (nova-one.org).

Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					
Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (C): Verification of Accuracy / Precision

PSA

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean	Expected Ranges
Level 1	2.68	1.57 – 3.78
Level 2	22.33	14.01 – >25.00
LOT #	6361A23002, 6362A23002	
Exp. Date	2025-11-30	

* Specifications vary by Lot#. For lot-specific values, please visit NovaOne's website (nova-one.org).

Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (C): Verification of Accuracy / Precision

Testosterone

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean	Expected Ranges
Level 1	393.4	221.0 – 565.9
Level 2	1307	819.1 – >1500
LOT #	6361A23002, 6362A23002	
Exp. Date	2025-11-30	

* Specifications vary by Lot#. For lot-specific values, please visit NovaOne’s website (nova-one.org).

Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					
Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (C): Verification of Accuracy / Precision

TSH

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean	Expected Ranges
Level 1	2.01	1.20 – 2.81
Level 2	17.82	11.18 – 24.46
LOT #	6361A23002, 6362A23002	
Exp. Date	2025-11-30	

* Specifications vary by Lot#. For lot-specific values, please visit NovaOne's website (nova-one.org).

Day 1 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					
Day 2 Perform 5 times Date ()	Level 1 Comment:					
Day 2 Perform 5 times Date ()	Level 2 Comment:					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (C): Verification of Accuracy / Precision

25-hydroxy vitamin D

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean	Expected Ranges
Level 1	40.89	21.53 – 60.26
Level 2	70.60	44.11 – >96.00
LOT #	6361A23002, 6362A23002	
Exp. Date	2025-11-30	

* Specifications vary by Lot#. For lot-specific values, please visit NovaOne’s website (nova-one.org).

<p align="center">Day 1 Perform 5 times</p> <p align="center">Date ()</p>	<p align="center">Level 1</p> <p>Comment:</p>					
<p align="center">Day 2 Perform 5 times</p> <p align="center">Date ()</p>	<p align="center">Level 2</p> <p>Comment:</p>					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**

Appendix (F): Temperature Log

Month _____ Year _____ Location _____

Day	Room Temp. 18 ~ 25 °C	Relative Humidity 20 ~ 80 %	Refrigerator 2 ~ 8 °C	Freezer < -20 °C	Comments/ Actions	Tech
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18						
19						
20						
21						
22						
23						
24						
25						
26						
27						
28						
29						
30						
31						

Appendix (H): External QC New Lot Number Verification

Perform this verification BEFORE you run out of the current control lot number. Before putting a new lot of control materials into use, verify that the new control performs as expected by parallel testing with the current lot number. All lot results should be within range each day. Run the new control five times, preferably on more than one day and using more than one operator. If the new control results are within range, you have verified the new control lot and can put it into use when the current QC lot is depleted. You will want to take action to determine why the control is not performing as expected when the new control results are outside of the range. Copy additional pages in case you needed.

	Level	Exp. Date (MM/DD/YY)	LOT #		Mean	Expected Range		
Current QC lot	Level 1	/ /				-		
	Level 2	/ /				-		
New QC lot	Level 1	/ /				-		
	Level 2	/ /				-		
Test item & Cartridge lot #		Date	Run	Level 1 result	Level 2 result	Tested by	Verified (Y/N)	Reviewed by
	Current QC lot	/ /	1					
		/ /	2					
		/ /	3					
		/ /	4					
		/ /	5					
	New QC lot	/ /	1					
		/ /	2					
		/ /	3					
		/ /	4					
		/ /	5					
	Current QC lot	/ /	1					
		/ /	2					
		/ /	3					
		/ /	4					
		/ /	5					
	New QC lot	/ /	1					
		/ /	2					
		/ /	3					
		/ /	4					
		/ /	5					

Appendix (I): Specimen Log

Date	Patient ID	Patient Name	Dr.	PSA Result	TSH Result	Free T4 Result	Testo. Result	Vit.D Result	Tech

Appendix (J): Internal Procedural Control Error Code Log

Every FREND™ Cartridge contains an Internal (built-in) Procedural Control feature. Each time a test is run in FREND™ System, the procedural zone is scanned and it validates the functional integrity of reagents. Detailed function of procedural control is explained in CLIA Compliance Manual. When FREND™ Cartridge’s procedural control is failed to show its validity, FREND™ System does NOT report the test result, instead, shows error code “EMTF-02, -04”.

Failure of internal procedural control might be caused from:

- improper sample volume,
- incompleteness of sample capillary flow,
- incompatibility caused from unknown interferences,
- damaged cartridge, and
- Mishandling of cartridges.

Technical support in NanoEntek provides full support on identifying the cause of an error when it frequently appears during patient sample test. You may record EMTF-02,-04 errors to monitor the frequency of QC failure during the test. Please follow a guide “How to Back-up & Send Error Data” to send an e-mail with error data for further investigation. To find out more regarding other error codes, refer to Maintenance section on the FREND™ System Manual.

Address	240 Bear Hill Road, STE 101. Waltham, MA 02451
E-mail	support@nanoentek.com
Tel.	781-472-2558 Toll Free: 888-988-0108

To check FREND™ System version, go to Setup → System View

FREND™ Serial Number	F10U		
Software Version		Firmware Version	

Test Date	Cartridge LOT#	Patient ID	Error Code	Repeat ¹ (Y / N)	QC Cartridge ² (Pass / Fail)	External QC run Level 1 ³ (In / Out)	External QC run Level 2 ³ (In / Out)

¹ Re-test the same sample one more on another cartridge to check if the same error occurs. (Yes = Same Error code occurred)

² Run QC Cartridge to check system integrity. Write “Pass” if QC is passed.

³ Run Level 1 and 2 of external QC solution from CLINIQA. If the results are within the range, write “in”
Check the record of QC if it passed or failed. If all results from external QC test are within the range, it is ‘Pass’

Please attach this Log with the error back-up data from FREND™ System when sending an e-mail.

Appendix (K): How to Back-up & Send Error Data

This is a guide to backup FREND™ System data to USB flash drive. Section A explains how to back up the data from FREND™ System and Section B will guide you how to compress folders for the attachment on your e-mail to technical support team.

You will need,


1. USB Flash drive for FREND™
2. Desktop or laptop which has a USB port
 - A. Windows / Mac OS (This example is prepared from Windows 10)
 - B. Support all OS versions
3. FREND™ System

Section A. Back-up Data from FREND™ System

1. Before start, back-up all your data in FREND™ USB flash drive into your computer to secure all the data.

**NOTE:* FREND™ System ONLY saves the latest sets of data.

2. Turn the FREND™ System on. Connect USB flash drive to the USB H on the back of the FREND™.



3. Click "**DATA**" in the menu and then click "**All**" button in the **Backup DATA** section.




**NOTE:* Do NOT click All in Delete DATA section. It will erase all data.

4. After data back-up is completed, press OK to exit and disconnect USB flash drive.


Appendix (K): How to Back-up & Send Error Data-continued

Section B. Send Backup Data to Technical Support

1. Connect USB flash drive to the computer. Find the "GIGASTICK (D:)" connected on your computer.




2. Open "GIGASTICK (D:)" and find two folders "CSV" and "Data"





* List of files and folders shown in this list may vary depending on your USB flash drive, but you can find CSV and Data folders easily if the backup was successful from Section A.

3. Select two folders and click "Share" tap on top of the window.



Appendix (K): How to Back-up & Send Error Data-continued

4. Find **"Zip"** button to compress two folders into one. New **"###.zip"** file will appear.



E-mail to:
support@nanoentek.com
For analysis of Errors

5. Attach this **"###.zip"** file to **support@nanoentek.com** with the Error Code Log records.

**NOTE:* The name of .zip file can vary depending on your computer.

Appendix (L): Incident Management Investigation Report Form

Attach copies of all original reports, amended reports, instrument printouts, and any other documentation as necessary.

Person Reporting Incident:	Reporting Date	Incident Date:
----------------------------	----------------	----------------

Information of Person Affected by Incident						
Affected Person:			DOB:		Patient Staff Other	
Ordering Physician:			Patient ID:			
Current Status:	Okay	Receiving Treatment	Hospitalized	Deceased	Unknown	
Notified	Name				Date Notified	
Physician						
Regulating Agency Ordering Physician						
Other Person						

Incident Description		
Error in phlebotomy or specimen collection procedure Misidentification of specimen Misidentification of patient Injury due to treatment	Misdiagnosis or injury based on test result Medication error based on test result Death based on test result Death Related to Treatment	Physical attack or abduction Failure in safety procedure Failure of instrument or methodology
Incident Description:		

True Cause Investigation Findings	
Factors	Underlying Factors
Human Equipment Controllable Environment Uncontrollable environment Other Factors_____	HR Issues Information Management Failure of Emergency responses Leadership Issues Other Factors_____
True Cause Investigation Findings:	

Appendix (L): Incident Management Investigation Report Form -continued

Patient Outcome	
Patient Outcome:	

Action Plan	
Opportunity for improvement	
Action to reduce reoccurrence	
Person Responsible for Implementation and Communication to Staff	
Implementation and communication Date	
Expected Results	

Action Plan Review	
Reviewer:	Review Date
Review Findings:	




Follow Up Actions		
No action required Tag Incident for Monitoring Future Evaluation of Corrective Action	Staff training and in-service competency assessment New policy or procedure Revision of policy or procedure	Cease Patient Testing Refer Patient Testing Resume Patient Testing

Reported By: _____ Reviewed By: _____

Follow Up Review By _____ Date Closed: __

Appendix (M): Laboratory Quality Incident & Correction Action

Attach copies of all original reports, amended reports, instrument printouts, and any other documentation as necessary.

Test System Affected:		Incident Date:
Analyte/Control Product:		Report Date:
Pre-Analytical	Analytical	Post-Analytical
		
Incident Description:		
Corrective Action Taken:		
Special Actions Taken:		
Physician Notified Lab Dir. / TC Notified	Specimen Rejected Specimen Recollected	Test Result Amended Test Canceled
Was there a procedural change as a result of this occurrence?		
Assessment and Action taken to Prevent Recurrence:		
Lab Director or Designee:		
Tag Incident for Monitoring	Future Evaluation of Corrective Action	Update SOP

Reported By: _____ Date: _____

Reviewed By: _____ Date: _____

Date Closed: _____

Appendix (N): Corrective Action Log

Serial# F10U _____ Year

Date	Problem Encountered	Corrective Action Taken	Tech

Reviewed by: _____ Date: _____

Appendix (O): Maintenance Log

Daily	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Wipe dust on FRENDO System																																
Power-on FRENDO System & View Screen																																
Record Temperature / Humidity																																
Monthly																																
USB Data Back-up																																
External QC																																
As needed																																
External QC New shipment, New lot																																
Internal procedural control Error Log																																

Reviewed by: _____, Date: