

Appendix (C): Verification of Accuracy / Precision

Vitamin D

Cartridge Lot#: _____ / Expiration Date: _____

External QC range		
	Mean(ng/mL)	Expected Ranges
Level 1		
Level 2		
LOT #		
Exp. Date		

* Specifications vary by Lot#. For lot-specific values, please visit CLINIQA's website (www.cliniqa.com), click the "Document Search" tab, enter your Lot number, and you will find the document containing the correct values.

Day 1 Perform 5 times Date ()	Level 1 Comment:					
	Level 2 Comment:					
Day 2 Perform 5 times Date ()	Level 1 Comment:					
	Level 2 Comment:					

Please, scan this file as PDF or Image and send to support@nanoentek.com or fax this document to **(781) 790-5649**