

## Appendix (C): Verification of Accuracy / Precision

# Free T4

Cartridge Lot#: \_\_\_\_\_ / Expiration Date: \_\_\_\_\_

External QC range		
	Mean(ng/dL)	Expected Ranges
<b>Level 1</b>		
<b>Level 2</b>		
<b>LOT #</b>		
<b>Exp. Date</b>		

\* Specifications vary by Lot#. For lot-specific values, please visit CLINIQA's website ([www.cliniqa.com](http://www.cliniqa.com)), click the "Document Search" tab, enter your Lot number, and you will find the document containing the correct values.

<b>Day 1</b> <b>Perform 5 times</b>  Date (            )	<b>Level 1</b>  Comment:					
<b>Day 2</b> <b>Perform 5 times</b>  Date (            )	<b>Level 2</b>  Comment:					
<b>Day 2</b> <b>Perform 5 times</b>  Date (            )	<b>Level 1</b>  Comment:					
<b>Day 2</b> <b>Perform 5 times</b>  Date (            )	<b>Level 2</b>  Comment:					

Please, scan this file as PDF or Image and send to [support@nanoentek.com](mailto:support@nanoentek.com) or fax this document to **(781) 790-5649**